

Disposable / Reusable Analog Forehead Thermometer Products

Instructions for Use

How to Read

1. Apply to center of forehead, just above the eyebrows and read temperature when yellow line stops moving, generally 15 seconds. **Read the indicator while it is still on the forehead.**
2. Center of yellow line indicates the correct temperature

Intended Use

- Measure human forehead surface temperature and monitor temperature trends
- After use, discard the disposable product or clean reusable product

Precautions

- Environmental temperatures might affect the performance of this product

Liquid Crystal thermometers are to be used on the forehead area only and estimate the core body temperature by measuring the temperature of the forehead. It is an ideal way to check for changes in temperature. A Liquid Crystal Forehead Thermometer is not as accurate for determining the actual core body temperature as a clinical (oral / rectal) thermometer. Instead it is designed to be an unbreakable, quick reading device for temperature screening between uses of a clinical thermometer.

Forehead skin temperatures are lower than oral temperatures.

Accurate to $\pm 1.0^{\circ}\text{F}$ / $\pm 0.6^{\circ}\text{C}$, the readings on a Liquid Crystal Forehead Thermometer have been adjusted to approximate oral temperature readings during normal to feverish conditions. If you think that your temperature is different from that shown by a Liquid Crystal Forehead Thermometer, you should check by using a clinical thermometer. If fever is present, consult your physician for appropriate treatment.

For External use only

This product contains no latex

Storage / Handling:

Products should be stored out of direct sunlight, dampness, or heat

Cleaning:

Wipe clean using soap and water or rubbing alcohol. Do not soak.



LCR Hallcrest, 1911 Pickwick Lane, Glenview IL 60026 USA



LCR Hallcrest Ltd. Dock Road, Connah's Quay, UK CH5 4DS

Notified Body: Det Norske Veritas **CE** 0434