


IIIB REUSABLE BABY BEAR BATH THERMOMETER CARD

Child Bath Time Safety is a Priority!



Safest temperature for bathing is 100°F.
American Burn Association

INSTRUCTIONS FOR USING BUBBLES BATH THERMOMETER:

- 1: There are 3 temperature indicators, Cool, Ok, Hot. Start with cool water and gradually add warm until "OK" is green. If no indicator lights, water is too hot or too cold.
- 2: Taking extreme care, thoroughly stir water to avoid potentially hot areas.
- 3: Immerse the bath thermometer completely for 15 seconds, ensuring that the water covers the black oval area of the thermometer. read the temperature of the thermometer while it is still in the water.
- 4: The correct temperature of the bath water will be indicated by the area on the thermometer scale that turns green.

This card is intended for use as a guide only and the manufacturers hold no responsibility for the actions of anyone based on the readings given

Markets:

- Hospitals
- New Parent Programs
- Daycares
- Family Wellness Promotions

Features:

- Reusable
- Helps Prevent Tap Water Burns in the Bath

Description:

IIIB Reusable Baby Bear Bath card with three event thermometer and custom imprint on front. Includes instructions for use on back. Start with cold water. Add hot, mixing well until the OK lights green. If no color appears the water is either too hot or too cold.

Specifications:

Overall Product Size: 2-1/8"W x 3-3/8"H x .020"
 Imprint Area Size: 1-3/4"W x 1/2"H
 Imprint on front.
 Temperature: Cool, OK, Hot
 91-94°F, 100-104°F, 113-117°F
 Weight: 8 lbs./M

Pricing:

Qty./Item	250	500	1000	2500	5000
IIIB	\$1.60	\$1.35	\$1.10	\$.92	\$.85

Price includes a one color imprint 5R

Other Charges:

Qty./Item	250	500	1000	2500	5000
Additional colors/imprint	\$.44	\$.30	\$.23	\$.20	\$.19

5T

Set up/Copy Change: \$31.25v

PMS Color Match Charge: \$25.00v

Production: 7-10 working days from art/order approval